



Green Klean® Chlorinated Disinfecting Tablets

A Better Choice for Disinfecting & Sanitizing With ALL Delivery Systems



All Microfiber
cloths



All Spray
Bottle



All Dry Wipes
Systems



All Microfiber Mop Heads



All Brands of
Electrostatic Sprayers



Buddy Jugs

HMIS Rating*
0-0-0

**6.55G Tablets –
Hospital grade disinfectant
120 tabs/6 bottles/case**

GK-CDT6.55

** In solution.*

It's easy to disinfect with Green Klean Chlorinated Disinfecting Tablets!

Tablets added to a quart of water meet the required ppm effective against a broad range of critical pathogens including C diff, C. auris, all COVID-19 viruses, HIV, Hepatitis A, Hepatitis B - Hepatitis C – TB – MRSA and more.

- No wait time for reentry with Electrostatic
- No Chemical Cartridge required with N95 Mask when used in electrostatic sprayers.
- No sticky residue left behind
- Broad Spectrum Disinfectant & Sanitizer
- No Pouring or Glugging
- Always Accurate - Safer
- 90% Less Storage Space
- Biodegradable
- Eliminates Quat Binding
- Replaces all Sanitizers, Disinfectants, Peroxides, and Bleach Products
- HMIS Rating in Solution 0-0-0
- Made in the USA

**Did You
Know?**

Our tablets are
Biodegradable
and the jars are
Recyclable!

Distributed by:



EPA # 71847-6-91038

It's easy to disinfect with Green Klean Chlorinated Disinfecting Tablets. Apply solution with a spray device, cloth, sponge, brush, or mop then allow to air dry, no rinse required.

Dwell Time	Pathogen	PPM	Dilution (Number of 6.55-gram Tablets)
1 Minute	Norovirus	2153	(1) tablet per quart
4 Minutes	C Diff	4306	(2) tablets per quart
2 Minutes	Candida Auris	4306	(2) tablets per quart
1 Minute	Aspergillus fumigatus	4306	(2) tablets per quart
4 Minutes	Candida albicans	4306	(2) tablets per quart
1 Minute	Coxsackievirus	4306	(2) tablets per quart
1 Minute	Hepatitis B and C	4306	(2) tablets per quart
2 Minutes	Hepatitis B and C	1076	(1) tablet per 2 quarts
10 Minutes	Hepatitis A	1076	(1) tablet per 2 quarts
1 Minute	Hepatitis A	4306	(2) tablets per quart
2 Minutes	HIV-1	1076	(1) tablet per 2 quarts
1 Minute	HIV-1	4306	(2) tablets per quart
4 Minutes	SARS-CoV-2 (COVID-19)	2153	(1) tablet per 1 quart
4 Minutes	Acinetobacter baumannii	4306	(2) tablets per quart
10 Minutes	VRE	1076	(1) tablet per 2 quarts
4 Minutes	VRE	4306	(2) tablets per quart
4 Minutes	CRKP	4306	(2) tablets per quart
10 Minutes	K pneumoniae	1076	(1) tablet per 2 quarts
10 Minutes	RSV	538	(1) tablet per gallon
10 Minutes	Rhinovirus (Common Cold)	1076	(1) tablet per 2 quarts
10 Minutes	Influenza (H1N1) Flu	1076	(1) tablet per 2 quarts
10 Minutes	Herpes	1076	(1) tablet per 2 quarts
10 Minutes	Poliovirus Type 1	1076	(1) tablet per 2 quarts

MARKETS



Education



**Sports
and More**



Senior Care



**Building
Service Cleaners**



Food Service



Hospitality



Bio-Security

****should not be used in a fogging device with less than 40 micron setting.***

Distributed by: